

Marc Chaigneau AFM/Raman Product Manager

High-speed Tip Enhanced Raman imaging

11-10-2016

Journée Thématique - Couplage Microscopies Optiques/Microscopie à Force Atomique Institut Pasteur de Lille

© 2015 HORIBA, Ltd. All rights reserved. Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

Bringing Scanning Probe Microscopies...

- SPM bring a lot of information on the physical characteristics of materials
 - Topography
 - Mechanical properties
 - Electrical and magnetic properties
- SPM is truly a nanoscale imaging technique...
 ...but it lacks *chemical* sensitivity



SS-DNA on HOPG functionalized with octadecylamine. 500nm frequency shift image.

HORIBA



AIST-NI

... and Raman together

Confocal Raman Microscopy is a very specific chemical imaging

- Precise structural information, wide areas of application
- Non-destructive technique, compatible with many environments
- A wide spectrum of available laser sources (from UV to IR : possibility of resonant Raman scattering)

Drawbacks

- Low cross-section (~ 10⁻³⁰)
- Limited spatial resolution



Graphene- HORIBA

- 156 x 180 = 28080 spectra (step = 0.5 μm)
- 2 min 08 (EMCCD, SWIFT; Acq. Time 2 ms + 1.5 ms)

© 2016 HORIBA, Ltd. All rights reserved.



Let's break the Rayleigh criterion!



AIST-NT



Nano Lightning Rod

plays the role of a Nano-Antenna

. Signal Enhancement

Near-Field Resolution



© 2016 HORIBA, Ltd. All rights reserved.



Conventional Raman VS TERS

Confocal Raman and TERS of the same area, graphene oxide and CNTs on Au



Confocal Raman 13 mW; integration 1 s





TERS 130 μW; integration 0.2 s

© 2016 HORIBA, Ltd. All rights reserved.







I. Bringing SPM and Raman together

II. Applications of TERS in Materials and Life Sciences

III. Stimulated Tip Enhanced Raman imaging

© 2016 HORIBA, Ltd. All rights reserved.



NanoRaman™





Bringing SPM and Raman together



XploRA Nano

© 2015 HORIBA, Ltd. All rights reserved.

Top and Side

Reflection Unit

NanoRaman: Optical configurations





- Top: co-localized measurements
- Side: optimized for TERS



AIST-NT

Horiba XploRA





Different approaches for the SPM



AFM: Atomic Force Microscope STM: Scanning Tunneling Microscope SF/NF: Shear-Force/ Normal Force Microscope → HORIBA TERS tips are available

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future



HORIBA





TRIOS Platform

Transmission and Reflection optical access



Versatile optical coupling platform

- Top-down, side, inverted access
- Ideal for custom optical components



© 2016 HORIBA, Ltd. All rights reserved.



TER signal comparison between Ag and Au tip

HORIBA

HORIBA



© 2015 HORIBA, Ltd. All rights reserved

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

Enhancement Factor of the Ag tips







© 2015 HORIBA, Ltd. All rights reserved.

HORIBA





- + Raman active layer: Ag with protective layers
- + Great enhancement factor **EF** = 10⁵⁻⁶
- + Easy-of-use thanks to the AFM regulation
- + Usable in top/bottom and side configurations
- + 9 tips out of 10 show the nanoresolution!

Excellent reproducibility of TERS maps





Tip #2







TERS: spatial resolution



Optical resolution capability: 8 nm Pixel step: 1.3 nm \rightarrow chemical sensitivity in both X and Y direction

Recorded @ SPIE San Diego, 2015

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific







I. Bringing SPM and Raman together

II. Applications of TERS in Materials and Life Sciences

III. Stimulated Tip Enhanced Raman imaging

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





<u>Why TERS</u>? → to give a fundamental insight into the **optical and electronic properties**



300 x 160 nm (100 x 60 points); Excitation 638 nm; 0.13 mW; integration 0.1 s

© 2016 HORIBA, Ltd. All rights reserved.



Raman at the

single tube

level

© 2016 HORIBA, Ltd. All rights reserved

Explore the future

1D materials: Nanotubes and Nanowires ALETINT

- Defect sites
- Discrimination between n- and p- doping
- Chirality variations (from the different radial breathing modes)

Quality, defect

density

- Transition from semiconducting to metallic behaviour
- Pressure induced frequency shift
- Decorated nano-objects
- Test sample for TERS





Impact on the

design of

devices

HORIBA

TERS on 2D materials





<u>Why TERS</u>? \rightarrow to characterize Graphene and 2D TMDCs in terms of size, shape, electronic properties, distribution of defects and contaminants



- a) 100 pixels per line TERS map of D-band intensity
- b) Topography image of the same flake
- c) representative TERS spectra
- d) distribution of the ratio of G to D band intensities

Increase of both the D and G TERS band intensities at creases and wrinkles

© 2016 HORIBA, Ltd. All rights reserved.





TERS Response of Graphene Oxide







2x0.75µm (125pixels/line) TERS map of GO on gold and current map recorded during the TERS map acquisition

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



TERS on 2D materials

<u>Why TERS</u>? → to characterize Graphene and 2D TMDCs in terms of size, shape, electronic properties, distribution of defects and contaminants

TEPL TERS TERS AFM image 200 nm AFM image 200 nm TERS

Y. Okuno, M. Chaigneau, HORIBA Scientific F. Fabbri, IMEM (unpublished).

- CVD grown MoS2 on <u>Si substrate</u>
- Ag tips
- Reflection configuration

Explore the future



TERS on 2D materials

<u>Why TERS</u>? → to characterize Graphene and 2D TMDCs in terms of size, shape, electronic properties, distribution of defects and contaminants



Kastl, Chen, Kuykendall, Darlington,

Borys, Krayev, Schuck, Aloni & Schwartzberg (LBL), Andrey Krayez (AIST-NT) (submitted).

© 2016 HORIBA, Ltd. All rights reserved.

Explore the future









- Identification of the defects
 (point defects, i.e. vacancy or dopants, or grain boundaries), edges
 Hydrogen-terminaned graphene (CVD)
 Local strain distribution, inhomogeneous and unintentional doping
- Decorated 2D materials



Graphene but also MOS2, WS2, WSe2, h-BN, GaSe

© 2016 HORIBA, Ltd. All rights reserved.



Just a Few Molecules

<u>Why TERS?</u> Single molecule spectroscopy is the **ultimate goal**!

Chemical mapping of C60 and C70 single molecules in gap mode AFM-TERS



© 2016 HORIBA, Ltd. All rights reserved.





- Single Molecule detection
- Organic molecules
- Polymer blends
- \rightarrow distinguish nanoscopic domains with distinct chemical signatures
- Photocatalysis, plasmon induced catalysis processes



PMMA, PEDOT/PSS, Polystyrene, Dye molecules, Thiols...

© 2016 HORIBA, Ltd. All rights reserved.





Semiconducting Nanostructures

<u>Why TERS</u>? Monitoring the local stress/strain in semiconducting nanostructures and the impact on the electronic properties





TERS spectra measured with the tip in tunneling position on the Si substrate (lower trace) and on the SiGe nano-stripe (upper trace).

© 2016 HORIBA, Ltd. All rights reserved.

Local stress measurement in nanostructures



© 2016 HORIBA, Ltd. All rights reserved.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

HORIBA

Life Sciences applications

Why TERS? To investigate isolated bio molecules (amino acid and nucleic acid)

Sequencing of Engineered Standards DNA on glass substrate • Ag tips Height(Dac) **Transmission configuration** Peak Filter 📕 🊖 💾 📖 🔻 5 nm Attr. Thymine (1395 cm-1) 0.2 10.4 um TERS spectral mapping of engineered DNA \rightarrow clear differentiation of spectral regions of - Whater hatte pattern and size consistent with the expected A/T and G/C homopolymeric blocks. 20 nm Attr. Cytosine (1280 cm-1)

courtesy of RMD, Noah Kolodjieski

© 2016 HORIBA, Ltd. All rights reserved.



Life Sciences applications

Correlation-based reconstruction of the original sequence



Sequence (image read bottom to top):

© 2016 HORIBA, Ltd. All rights reserved





- Bio-molecules
- Nucleic acids \rightarrow on the road to DNA sequencing.
- Amino acids (investigation of protein structures).
- Micro-organisms:
- Pathogens (bacteria, viruses) .
- Lipid membranes (cell surface).

Nanometre length-scale resolution of TERS Label-free measurement + fingerprint specificity Identification, detection, molecular dynamics HORIBA





I. Bringing SPM and Raman together

II. Applications of TERS in Materials and Life Sciences

III. Stimulated Tip Enhanced Raman imaging

© 2016 HORIBA, Ltd. All rights reserved.



Spontaneous VS stimulated TERS



- Enhancement factor = 10⁵⁻⁶
- 100 1000 cts/sec
- S/N = 10 50



expected Gain = 10⁹!





$$G = \frac{S_{STIM}}{S_{SPON}} = [32 \ \rho^3 c^2 / W_s^2] F(W_s)$$





Stimulated TERS setup



TERS spectrum and 'tuned - detuned' STERS of 1141 cm⁻¹ vibration of azobenzene

 $\ensuremath{\textcircled{}}$ 2016 HORIBA, Ltd. All rights reserved.



Mapping of patterned azobenzene thiol SAM



simultaneously recorded STM and STERS images *continuous* (raster scan) mapping: *14 min vs 2 h (theoretically) in spontaneous TERS!*

© 2016 HORIBA, Ltd. All rights reserved.

K. H. Wickramasinghe, M. Chaigneau, et al ACS Nano (2014)

HORIBA

HORIBA

STERS: spatial resolution

Azobenzene thiol SAM maps of smaller area (180 nm x 180 nm)



Simultaneous Images @ 0.5 Hz/lines

Optical resolution : < 6 nm

Stimulated gain : 109

 $S/N = 2.10^3$

K. H. Wickramasinghe, M. Chaigneau, et al ACS Nano (2014)

© 2016 HORIBA, Ltd. All rights reserved. Explore the future HORIBA

AJST-NT

CONCLUSIONS

TERS <u>*IS READY*</u> for real-life analytical application on number of scientifically/industrially interesting samples.

- Individual components in complex samples of carbon-based materials can be conventionally imaged <u>with about 10 nm</u> resolution and unambiguously identified based on their Raman spectra.
- Without detailed NANOSCALE characterization of molecules, 1-d and 2D materials their use in next generation devices is pretty much impossible. TERS can provide reliable information on peculiarities of local structure with resolution below 10 nm.

The field of TERS is now rapidly advancing past technique development toward answering important scientific questions in materials science, physics, surface spectroscopy and chemistry.

© 2016 HORIBA, Ltd. All rights reserved.



CONCLUSIONS

new trends

Non linear TERS: Stimulated TERS demonstrated

- Fast mapping capability
- Gain in SNR (with respect to TERS)
- New SPM chemical mode \rightarrow potential applications

UHV TERS

- **TERS in liquids**
- **Polarized Tip Enhanced Raman Scattering**
- **ElectroChemical Tip Enhanced Raman Spectroscopy**

New hybrid cantilever-based TERS tips

.....

© 2016 HORIBA, Ltd. All rights reserved.



Acknowledgement

Andrey Krayez Alexei Temiryazev Sergey Bashkirov Victor Baukov **Dmitry Evplov Vladimir Zhizhimontov Vasily Gavrilyuk Eddy Robinson Sergey Saunin Emmanuel Leroy Joachim Schreiber** Patrick Hsia **Damien Eschimese Ophélie Lancry Jérémy Brites Yoshito Okuno** Kumar Wickramasinghe Razvigor Ossikovski Ryohei Yasukuni





Thank you very much for your attention.

marc.chaigneau@horiba.com

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

